

Wavelength Range:	190 - 1100 nm 800 - 1800 nm 1.5 - 3.5 $\mu\text{m}$
Resolution:	2 $\mu\text{m}$
Smallest Beam:	100 $\mu\text{m}$
Scanned Area:	5 x 7 mm - Si 3 mm - Ge 3 mm - InGAs
With 2D Stage:	23 x 43 mm

# BeamScope-P8

## Typical Applications

- ◇ Laser/Diode Laser characterization
- ◇ Laser assembly development, alignment, characterization, production test & QA.
- ◇ Lasers And Laser Assemblies for
  - Disk/Wafer Characterization
  - Laser Printing/Marking
  - Medical Lasers
  - Bar Code Scanners ... etc.

## Real Time:




- ◇ X-Y Profile Measurement
- ◇ Angular Divergence
- ◇ Ellipticity, Centroid, Gaussian fit
- ◇ Relative Power

## Features:

- ◇ Beam dimensions 100  $\mu\text{m}$  to 45 mm
- ◇ Resolution 2  $\mu\text{m}$  or 0.5 %
- ◇ 190 nm to 3.5  $\mu\text{m}$  options
- ◇  $M^2$  measurement accessory
- ◇ ISO 11146 Compliant
- ◇ Narrow probe for confined spaces
- ◇ Front mounted apertures
- ◇ Wide dynamic range
- ◇ Powerful, intuitive software
- ◇ Upgrade option, BeamScope-P7 to USB 2.0

## Accessories & Options

- ◇  $M^2$  Measurement Accessory-USB 2.0
- ◇ UV - NIR Options 190 to 3.5  $\mu\text{m}$
- ◇ 2-D Scan Stage for 23 x 45 mm profiling



BeamScope-P8 system is comprised of: a compact head, interface box, a 3 m (10 ft.) USB 2.0 cable, user manual, and Software for Windows XP, Vista or Windows7. (Windows and Vista are registered trademarks of Microsoft Corp.)



**Principle of Operation:** A linear scanning probe carries either a single pinhole, a single slit, or orthogonal X-Y slits. This linear scan satisfies the strict requirements of the ISO 11146 laser profiling standard\*. Light passing through the slits falls onto a Silicon (190 to 1150 nm) or Germanium (800 to 1800 nm) detector. [\* Until the introduction of the DataRay BeamScope Beam Profiler, no commercially available slit scan or knife-edge scan beam profiler met the ISO 11146 Standard. The Standard requires that the scan be performed in a *plane* orthogonal to the propagation axis. Drum style scanners cannot meet the Standard. DataRay's unique linear scan probe is designed to fully comply with the Standard.]

## Acquire Beam Profiles In Constricted Areas

BeamScope-P8 has made the measurement of once inaccessible beam profiles not merely possible, but simple. The unique probe-style scan head easily peers into confined axial gaps between lens, mirrors, and filters. Its ability to probe along-axis spaces as narrow as 12mm creates a whole new world of applications.

## No Beam Distortion From Optics Or Filters

There's no distortion of the beam due to ancillary optics or filters because the BeamScope-P8 doesn't need them when analyzing most type of lasers. The AUTO GAIN feature can continually adjust the detector amplifier gain to ensure full use of the 55 dB (300,000:1) gain range. Spot dimensions from 3  $\mu\text{m}$  to 23 mm can be measured from a single scan head. Scan beam areas up to 23 x 45 mm with the new 2-D stage accessory.

## Front Mounted Apertures

Front mounted apertures enable you to see precisely where the beam is being measured. Rapidly diverging and fast focusing beams are simple to capture if you can get close to your source. Now it's easier than ever to measure laser diode arrays, micro-lensed sources, broad stripe lasers, etc. Available apertures can accommodate power densities up to 100 W/mm<sup>2</sup> onto small pinhole apertures. (Max. Total power = 0.5 W) Change apertures in minutes from slits to pinholes. This makes the BeamScope-P8 an unbeatable value in beam analyzers.

## Notebook PC Portability with USB 2.0


BeamScope plugs interface to a USB 2.0 port on your notebook PC to give a portable unit with a small footprint.

## Perfect For R&D, QA & Production

R&D users will appreciate the comprehensive range of analytical functions. QA & Production engineers will appreciate the ability to save test configurations as JOB files, and to indicate **Pass/Fail** on-screen.


## Unparalleled M<sup>2</sup> Measurements

The BeamScope™ P8 optional M2DU-P8 accessory is unlike any other on the market. There are no complicated adjustments, yet the user can achieve highly repeatable measurements. The software automatically takes more frequent measurements in the waist area in order to accurately determine the true waist diameter.



## 23 x 45 mm 2-D Scan Stage

DataRay's latest beam profiling innovation offers *extraordinary* 0.2% (512 x 512 pixels) resolution, down to 5 x 5  $\mu\text{m}$  (HxV), over beam areas up to 23 x 45 mm. 2-D scan results display in the integral imaging software for area image analysis.



<b>Measurable Sources</b>	CW or Pulsed sources $\geq$ 5kHz Pulse Rep Rate @ 5% duty factor. Higher PRR is better.
<b>Measured Beam Powers</b>	See the graph in the Notes, below. E.g. 6 $\mu$ W to 3 W, for a 1 mm diameter ( $1/e^2$ ) Gaussian beam @ 633 nm, 5 $\mu$ m slit.
<b>Optical Dynamic Range</b>	55 dB (= 300,000:1) [75 dB with Neutral Density 2.0 film]
<b>Shape of Maximum Scanned Area</b>	<b>Important:</b> For accurate measurements, beam width should be $\leq 0.5 \times$ Scan Dimensions With /EPH extended probe head, dimension 23 mm below, becomes 35 mm.
<b>Pinholes (PA series)</b>	<b>Shape</b> <u>Cross Scan x Scanned Length</u>
<b>Single Slits (SS series)</b>	Line Scan Pinhole diameter x 23 mm
<b>X-Y Slits (XY series)</b>	Rectangle 7* x 23 mm, (* 5 for Ge, 3 for InAs)
<b>2-D Stage (M2B)</b>	Trapezoid 5* x 15/5 mm, (* 3 for Ge, 2 for InAs, 3.5 x 13.5/6.5 for XYPI5)
	Rectangle 45 x 23 mm scanned area image. Scans a pinhole over this area.
<b>Measured Beam Diameters/Widths</b>	100 $\mu$ m to $\sim$ 25 mm. Defined as the $1/e^2$ diameter, = 13.5% of peak for Gaussian beams)
<b>Measurement Resolution</b>	2 $\mu$ m, or 0.5 % of the measured beam diameter, whichever is greater
<b>Measurement Accuracy</b>	$\pm 2 \mu$ m $\pm$ 2% of measured beam diameter
<b>Measured Beam Profiles</b>	X & Y Linear & logarithmic profile display modes
<b>Measured Profile Parameters</b>	Gaussian beam diameter Gaussian fit Second Moment beam diameter Knife-Edge beam diameter Centroid position, relative and absolute Ellipticity + Orientation of Major Axis Beam Wander display
<b>Displayed Profiles</b>	X only, Y only, X & Y 2-D plot (10,16 or 256 colors) reconstruction see Note 1 3-D plot (10,16 or 256 colors) reconstruction see Note 1
<b>Update Rate</b>	1 to 2 Hz. Depends upon the PC Processor Speed, Scanned Profile & Selected Options
<b>Data Analysis</b>	
<b>Pass/Fail</b>	On all measured parameters, on-screen, in selectable <b>Pass/Fail</b> colors
<b>Averaging</b>	Beam Diameter Running Average and Accumulating Average options
<b>Standard Deviation</b>	On Accumulation Averaging Screen
<b>Power Measurement</b>	Units of mW, dBm, dB, % or user entered (relative to a reference measurement provided by the user.)
<b>Source to Slit Distance</b>	1.0 mm minimum
<b>Aperture sizes</b>	<b>Important:</b> See Scanned Area (above) for measurable beam dimensions see Note 2
Slits	2.5, 5, 10, 25 and 100 $\mu$ m wide 7 mm long (Planar version of 5 $\mu$ m slits are 5 mm long)
Pinholes	5, 10, 25 and 50 $\mu$ m diameter (Larger or smaller pinholes to special order) see Note 2
<b>Wavelength Range</b>	
Silicon Detector	190 to 1150 nm
Germanium Detector	800 to 1800 nm
InAs	1500 to 3.5 $\mu$ m
<b>Mounting</b>	1/4-20 & M6 threaded mounting holes
<b>Temp. Range (inc. Accessories)</b>	
Operating	10° to 35° C
Storage	5° to 45° C
<b>Minimum PC Requirements PC or Intel-Mac</b>	USB2.0 port, Windows XP, Vista or Windows 7; 1024 MB RAM; 10 MB Hard Drive space; 1024 x 768 monitor. Windows and Vista are trademarks of Microsoft Corp.



## Notes:

1. The 2-D and 3-D profiles are 'reconstructed' from the X-Y scan, making the assumption that the measured X beam profile is the same for all values of Y, and that the measured Y beam profile is the same for all values of X.
2. In Single Slit or Pinhole mode, the slit/pinhole width should be  $\leq 1/3^{\text{rd}}$  of the diameter of the beam under measurement. For X-Y slit pairs inclined at  $\pm 45^\circ$ , the ratio is approximately  $1/5^{\text{th}}$ .
3. In Knife-Edge mode, the slit width should be  $\geq 3x$  the beam diameter. For X-Y slit pairs inclined at  $\pm 45^\circ$ , the ratio is  $\geq 4.3$  times the beam diameter.

**Power Limit.** The graph allows you to simply determine the approximate maximum optical power that BeamScope can measure without additional attenuation. The limit is a detector current limit.



### PxS is the Power Limit in Watts.

To calculate the power limit at the laser wavelength:

The graph shows the approximate saturation beam power in mW versus wavelength for a 100  $\mu\text{m}$  diameter beam ( $1/e^2$ ) & a 5  $\mu\text{m}$  slit.

- 1) Determine the plotted value,  $P_S$  mW, for your wavelength. (e.g. 70 mw at 633 nm)
- 2) For a different slit or beam, multiply  $P_S$  by:  
 $0.05 \cdot (\text{beam diam., } \mu\text{m}) / (\text{slit width, } \mu\text{m})$
- 3) For a pinhole multiply  $P_S$  by:  
 $0.05 \times (\text{beam diam., } \mu\text{m})^2 / (\text{pinhole diam. } \mu\text{m})^2$

Approximate saturation beam power  
100  $\mu\text{m}$  beam diameter, 5  $\mu\text{m}$  slit.



## Part Numbers

## BeamScope Products & Accessories

### System Examples

[Includes Software, 3 m long USB 2.0 cable, User Manual]

<b>BS8-XY5</b>	Silicon Detector system with 5 $\mu\text{m}$ X-Y Slit
<b>BS8-XY2.5</b>	Silicon Detector system with 2.5 $\mu\text{m}$ X-Y Slits
<b>BS8G-XY5</b>	Germanium Detector system with 5 $\mu\text{m}$ X-Y Slits
<b>BS8G-PA10</b>	Germanium Detector system with 10 $\mu\text{m}$ Pinhole
<b>BS8-InAs SXX</b>	Indium Arsenide Detector system with Single slit (pinholes are also available)
<b>/EPH</b>	Suffix for +10 mm Extended probe head, for probing deeper recesses P8 USB 2.0 interface

### Accessories

<b>PA5</b>	Pinhole 5 $\mu\text{m}$ diameter
<b>PA10</b>	Pinhole 10 $\mu\text{m}$ diameter
<b>PA25</b>	Pinhole 25 $\mu\text{m}$ diameter
<b>PA50</b>	Pinhole 50 $\mu\text{m}$ diameter
<b>PA100</b>	Pinhole 100 $\mu\text{m}$ diameter
<b>SS2.5</b>	Single Slit 2.5 $\mu\text{m}$ wide x 3 mm long
<b>SS5</b>	Single Slit 5 $\mu\text{m}$ wide x 7 mm long
<b>SS10</b>	Single Slit 10 $\mu\text{m}$ wide x 7 mm long
<b>SS25</b>	Single Slit 25 $\mu\text{m}$ wide x 7 mm long
<b>SS50</b>	Single Slit 50 $\mu\text{m}$ wide x 7 mm long
<b>SS100</b>	Single Slit 100 $\mu\text{m}$ wide x 7 mm long
<b>XY2.5</b>	X-Y Slit 2.5 $\mu\text{m}$ wide x 3 mm long @ $\pm 45^\circ$ Planarity $\pm \sim 40 \mu\text{m}$
<b>XY5</b>	X-Y Slit 5 $\mu\text{m}$ wide x 7 mm long @ $\pm 45^\circ$ " "
<b>XY10</b>	X-Y Slit 10 $\mu\text{m}$ wide " " "
<b>XY25</b>	X-Y Slit 25 $\mu\text{m}$ wide " " "
<b>XY50</b>	X-Y Slit 50 $\mu\text{m}$ wide " " "
<b>XY100</b>	X-Y Slit 100 $\mu\text{m}$ wide " " "
<b>XYPI5</b>	Dual planar X-Y Slit 5 $\mu\text{m}$ wide x 5 mm long, planarity set to $\pm 4 \mu\text{m}$ (for tightly focused beams).
<b>M2DU-BS</b>	M <sup>2</sup> Accessory, 44 mm scan / 2.5 $\mu\text{m}$ steps, + 3 m long USB cable + M <sup>2</sup>
<b>LNZ-UV-XXX</b>	Lens 185-450 nm, focal lengths from 50 mm to 1000 mm, clear apertures 22 and 47 mm
<b>LNZ-VIS-XXX</b>	Lens 400-700 nm, focal lengths from 50 mm to 1000 mm, clear apertures 22 and 47 mm
<b>LNZ-NIR-XXX</b>	Lens 630-1100 nm, focal lengths from 50 mm to 1000 mm, clear apertures 22 and 47 mm
<b>LNZ-TEL-XXX</b>	Lens 1030- 1800 nm, focal lengths from 50 mm to 1000 mm, clear apertures 22 and 47 mm
	Additional lenses/ focal lengths are available. Please contact us for custom configurations
<b>M2DU-2D</b>	2D Scanning Stage 44 mm travel, 2.5 $\mu\text{m}$ steps. Beam size to 23 x 44 mm with <b>M2DU-AP</b> plate
<b>M2DU-AP</b>	Additional adapter plate to connect BeamScope to M2DU-2D Scanning stage
<b>P7-P7U-UPGRD</b>	Upgrade a PCI BeamScope-P7 to USB 2.0 interfaced BeamScope-P7U with rebuilt head, calibration certificate, Interface box, USB 2.0 cables + Universal wall socket power supply and 3 year warranty.
<b>P7-USB2-IF</b>	Add USB2.0 interfacing to a BeamScope-P7.. [No head rebuild, no calibration, no warranty extension.]

## M2 Accessory / 2D Stage Specifications

<b>Standard Lens Assembly</b> (Detached for 2-D Stage operation)	Coated achromatic lens, $\lambda/4$ @ 550 nm ( see below for other wavelengths) Focal Lengths 50 mm to 500 mm, wavelengths UV to 16 microns Clear Aperture 22 and 47 mm
<b>Stage Travel</b>	44 mm
<b>Step Size</b>	2.5 $\mu\text{m}$
<b>Dimensions M2DU</b>	Across axis width x Height with BeamScope x Along axis depth 90 x 100 x 200 mm (5.5 x 6.5 x 8 inches)
<b>Weight M2DU</b>	0.8 kg (1.8 lb)
<b>M2DU + BeamScope-P8</b>	1.5 kg (3.3 lb)

## Other DataRay Beam Profiling Instruments

<b>BeamMap2</b>	Real Time M-Squared Multi-plane profiler 0.1 micron resolution on CW lasers. XYZ profiles, XYZ focus, Centroid, Alignment, Divergence, M <sup>2</sup> , Visible to Telecom wavelengths. Port-powered USB2.0
<b>Beam'R2</b>	0.1 micron resolution on CW lasers, 0.5 micron to 4 mm beam dimensions. Port-powered USB2.0
<b>WinCamD</b>	High performance, high resolution CCD and CMOS imaging systems. Port-powered USB2.0
<b>WinCamD-XHR</b>	New High resolution CCD array with 3.2 $\mu\text{m}$ pixel and 2048 x 1536 array Port Powered USB2.0
<b>WinCamD-FIR</b>	Wavelength Range 2-16 $\mu\text{m}$ , 25 $\mu\text{m}$ pixel pitch, 320 x 240 pixels, 7.5 x 6.6 mm active area
<b>BladeCamXR</b>	The Smallest Beam Profiler in the known Universe - 1/2" and 1/1.8" formats Pixel size 4.4 and 5.2 $\mu\text{m}$

