



**DataRay Inc.**  
**Calibration RMA Form**

Beam Profiling ... Engineered as a System  
 Delivered as a Solution!

www.dataray.com  
 Ver. 0710A

**Calibration:** Dataray offers annual recalibration on all BeamMap series and BeamScope P5/P7 products. Contact the factory or your local distributor for current pricing.

**BeamScope-P7:** The linear probe mechanism does wear with intensive use, therefore annual recalibration is recommended for instruments used regularly in production test.

**BeamMap series:** Annual recalibration is offered but is not considered necessary unless the head has been subject to physical abuse. The only moving part is the puck on a stepper motor. The ball bearing motor is good for 5 years continuous operation, and to date we have yet to see a system where motor wear is an issue. The calibration is held in an EEPROM in the head.

**WinCamD series:** No calibration is required. [Calibration is intrinsic to the software and the chip photolithography.]

**RETURN TO:** Important. TO SAVE TIME, return products to the CORRECT address:

**BeamScope:** DataRay Bella Vista Office, 14262 Troy Hills Drive, Bella Vista, CA 96008-0247, USA 530-472-1872

**BeamMap:** DataRay Colorado Office, 3840 Barr Court, Boulder, CO 80305, USA 303-543-8235

Fill out all relevant items marked with an asterisk (\*)

\* **Calibration RMA #:** \_\_\_\_\_ Obtain from [support@dataray.com](mailto:support@dataray.com) or call 303-543-8235

\* **Contact Name(s)** \_\_\_\_\_ [include a technical contact]

\* **Organization & Dept.** \_\_\_\_\_

\* **Ship-to Address** \_\_\_\_\_

\* **Phone #(s)** \_\_\_\_\_ **Fax:** \_\_\_\_\_ **Email:** \_\_\_\_\_

\* **Returned Items (X)** ALWAYS include head, cable & PC card.

	Head (Serial #)	PC Card	Cable(s), Quantity	Obvious Damage?
Example line:	X 1234	X	X 1	No
BeamMap-series	_____	_____	_____	_____
BeamScope	_____	_____	_____	_____

\* **Calibration Price:** \_\_\_\_\_ As quoted to by your Distributor or DataRay.

Name & Position: \_\_\_\_\_ PO No.: \_\_\_\_\_ **MUST PROVIDE**

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

Credit Card: \_\_\_\_\_ Acct.#: \_\_\_\_\_ Expires: \_\_\_\_\_ Name on Card: \_\_\_\_\_

\* **Return Shipment:** You **must** provide a FedEx (preferred) or UPS account number for return shipping.

FedEx \_\_\_ UPS \_\_\_ Other (specify) \_\_\_\_\_ Account No.: \_\_\_\_\_

Requested delivery: \_\_\_ Priority overnight \_\_\_ Standard overnight \_\_\_ 2-day \_\_\_ 2-day International

---

**DataRay use only.** Copy to customer, original to customer file. Calibrated by: \_\_\_\_\_ Date: \_\_\_\_\_ Charge: \$ \_\_\_\_\_

Comments: \_\_\_\_\_

---